


ZYSENSE MGA 800 ANALYZER

# ZYSENSE MULTI-GAS ANALYZER MGA -800



ZYSENSE MGA 800 - TABLET SCREEN



REAL TIME MEASUREMENT OF 5 GASES CONCENTRATION, FLOW & PRESSURE

## **ZYSENSE MULTI-GAS ANALYZER MGA- 800**

Zysense Multi-Gas Analyzer was developed at the request of several researchers in the Nitric Oxide research community to measure other two cell signaling molecules of interest, Hydrogen Sulfide, Carbon Monoxide. In addition to the above several other medical, bio pharmaceutical, and doctors in clinical research requested several others gases like Ammonia, Hydrogen, Phosphine, Nitrogen Dioxide, Formaldehyde and others. Since most of the gases interact and interfere with the measurements at low concentrations, it was critical to develop sensors that were accurate and were not influenced by the presence of other gases.

The medical research community also requested a versatile, robust equipment where they could change sensors and measure multiple parameters critical to a specific disease or measure certain exhaled gas markers that provide some “non invasive” insights into the health of the subject.

### **The unique features of Zysense MGA -800 are as follows:**

1. Unique custom developed sensors designed for clinical research applications ( Nitric Oxide, Nitrogen Dioxide, Carbon Monoxide, Hydrogen Sulfide, Ammonia, Hydrogen, Oxygen, Phosphne, and Formaldehyde).
2. Sensors are uniquely designed to minimize or eliminate interaction between the various gases that are being measured.
3. All gases can be measured simultaneously using the same instrument
4. Can be used to measure gases in a gas stream or exhaled breath, or from nasal exhalation from a subject.

### **Applications:**

5. Can also be used in neonatal incubator ( isolette) to measure exhaled gases.
6. Can be used to measure breath exhaled gases in sleep studies, and stress tests on a treadmill.
7. Can be used to measure the Nitric Oxide, Nitrogen Di Oxide, Oxygen levels while dosing nitric oxide to clinical subjects

### **Major Research currently underway that utilize Multi-Gas Analyzer Zysense MGA-800**

1. Exhaled breath for the measurement of microbiome in lower airway in lungs - looking for anaerobic conditions
2. Exhaled breath gas measurement during Pulmonary Function Tests ( Planned Clinical Trial in 2Q 2020)
3. Exhaled breath gas measurement for Cystic Fibrosis in children ( Planned trial)

**Potential advantages and limitations of exhaled breath analysis:**

**Advantages**

Non-invasive and non-intrusive  
Allows for repeated measurements  
Can be inexpensive  
Potential for portability  
Potential for real-time results

**Limitations**

Confounders such as diet and environment  
Lack of standardization  
Poor reliability of prototypes  
Storage  
Physician acceptance

**Mitigation Steps**

Managed DOE  
Standardizing methods  
Reliable Sensors ( MGA 800)  
Extremely rugged ( MGA-800)  
Multiple Clinical Trials  
with( MGA 800) by doctors

Personalized medicine (“breath-print”)

None

Future MGA 800 use in  
clinics





## HYDROCARBON SENSOR


Hydrocarbon vapor sensor is an electrochemical device used for the detection of a variety of hydrocarbon vapors in ambient air. This sensor is used for measurement of ethylene oxide, formaldehyde, alcohols, acetylene, ethylene, and others. The cross-sensitivity table given below lists the hydrocarbons that may be detected with this sensor. It is designed to be used in conjunction with Zysense MGA 800 Multi-Gas Analyzer. The HC sensor is a 3-electrode electrochemical cell and operates by generating a small electrical current proportional to the partial pressure of hydrocarbon vapor in the surrounding air. The current results from the direct oxidation of the hydrocarbon on the measuring electrode to form carbon dioxide. Oxidation of formaldehyde is shown in the formula below.




### Hydrocarbon Sensor - Cross Sensitivity Data

Gas	Symbol	Response to 1 PPM
Ethylene Oxide	C <sub>2</sub> H <sub>4</sub> O	1.0
Formaldehyde	HCOH	1.0
Methanol	CH <sub>3</sub> OH	1.2
Ethanol	C <sub>2</sub> H <sub>5</sub> OH	0.7
Isopropyl Alcohol	C <sub>3</sub> H <sub>7</sub> OH	0.25
Acetylene	C <sub>2</sub> H <sub>2</sub>	1.2
Ethylene	C <sub>2</sub> H <sub>4</sub>	1.0
Vinyl Chloride	C <sub>2</sub> H <sub>3</sub> Cl	0.6
Propylene Oxide	C <sub>3</sub> H <sub>6</sub> O	0.7
Propylene	C <sub>3</sub> H <sub>6</sub>	0.6
Methyl Methacrylate	C <sub>5</sub> H <sub>8</sub> O <sub>2</sub>	0.25
Acetaldehyde	CH <sub>3</sub> CHO	0.25
Styrene	C <sub>8</sub> H <sub>8</sub>	0.25
Carbon Monoxide	CO	0.5
Nitrogen Dioxide	NO <sub>2</sub>	0.1
Sulfur Dioxide	SO <sub>2</sub>	0.4
Methyl Mercaptan	CH <sub>3</sub> SH	1.0
Hydrogen	H <sub>2</sub>	0.05
Nitric Oxide	NO	0.8

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for an HC sensor exposed to Formaldehyde. Response to other hydrocarbons may be somewhat different. Note that response times can become significantly slower at temperatures below -20°C.



## CARBON MONOXIDE SENSOR

CO sensors exhibit response to certain other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any of these potential interferences are present and might present interference issues.

Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.

Carbon Monoxide sensor is an electrochemical device used for the detection of CO gas leaks in ambient air. It is designed to be used in conjunction with ZYSENSE MGA 800. Carbon Monoxide sensors operate by generating a small electrical current proportional to the partial pressure of CO gas in the surrounding air. The current is the result of the electrochemical oxidation of CO on a catalytic electrode as shown below. Carbon monoxide sensors are 3-electrode sensors and **require oxygen levels above 5% to operate** properly. They **may not be used** in **oxygen free** applications.




### Carbon Monoxide sensor cross sensitivity data


Gas	Symbol	Response to 1 PPM
Hydrogen	H <sub>2</sub>	0.1
Ethylene	C <sub>2</sub> H <sub>4</sub>	0.1
Hydrogen Cyanide	HCN	0.1
Acetylene	C <sub>2</sub> H <sub>2</sub>	0.1
Nitric Oxide	NO	0.1
Carbon Dioxide	CO <sub>2</sub>	None
Acid Gases	HX	None
Hydrogen Sulfide	H <sub>2</sub> S	None
Nitrogen Dioxide	NO <sub>2</sub>	None
Sulfur Dioxide	SO <sub>2</sub>	None
Phosgene	COCl <sub>2</sub>	None
Chlorine (Oxidants)	Cl <sub>2</sub>	None



Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for a carbon monoxide sensor. Note that this response time can become significantly slower at temperatures below -20°C.



## NITRIC OXIDE Sensor

Nitric Oxide sensor is an electrochemical device used for the detection of NO gas leaks in ambient air. It is designed to be used in conjunction with Zysense Multi-Gas Analyzer MGA - 800.

Nitric Oxide sensors operate by generating a small electrical current proportional to the partial pressure of NO gas in the surrounding air. The current is the result of the electrochemical oxidation of NO on a catalytic electrode as shown below. Nitric Oxide sensors are 3-electrode sensors and **require oxygen levels above 5% to operate properly**. They **may not be used in oxygen free applications**.




NO sensors exhibit response to certain other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any of these potential interferences are present and might present interference issues. Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.


### Nitric Oxide Sensor - Cross Sensitivity data:

Gas	Symbol	Response to 1 PPM
Hydrogen	H <sub>2</sub>	None
Ethylene	C <sub>2</sub> H <sub>4</sub>	None
Hydrogen Cyanide	HCN	None
Acetylene	C <sub>2</sub> H <sub>2</sub>	0.01
Phosphine	PH <sub>3</sub>	None
Carbon Dioxide	CO <sub>2</sub>	None
Acid Gases	HX	None
Hydrogen Sulfide	H <sub>2</sub> S	None
Nitrogen Dioxide	NO <sub>2</sub>	None
Sulfur Dioxide	SO <sub>2</sub>	None
Phosgene	COCl <sub>2</sub>	None
Chlorine (Oxidants)	Cl <sub>2</sub>	None

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for a nitric oxide sensor. Note that this response time can become significantly slower at temperatures below -20°C.



## NITROGEN DIOXIDE ( NO2) SENSOR

Nitrogen Dioxide sensor is an electrochemical device used for the detection of NO<sub>2</sub> gas leaks in ambient air. It is designed to be used in conjunction with Zysense Multi-Gas Analyzer MGA -800. .

Nitrogen Dioxide sensors operate by generating a small electrical current proportional to the partial pressure of NO<sub>2</sub> gas in the surrounding air. The current is the result of the electrochemical reduction of NO<sub>2</sub> on a catalytic electrode as shown below. Nitrogen Dioxide are 3-electrode sensors and require **oxygen levels above 5% to operate properly**. They **may not be used in oxygen free applications**.




H10-26 NO<sub>2</sub> sensors exhibit response to certain other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any of these potential interferences are present and might present interference issues. Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.


### Nitrogen Dioxide Sensor - Cross Sensitivity data

Gas	Symbol	Response to 1 PPM
Hydrogen	H <sub>2</sub>	-0.1
Ethylene	C <sub>2</sub> H <sub>4</sub>	-0.01
Hydrogen Cyanide	HCN	-0.02
Acetylene	C <sub>2</sub> H <sub>2</sub>	-0.02
Phosphine	PH <sub>3</sub>	-3.5
Carbon Dioxide	CO <sub>2</sub>	None
Hydrogen Chloride	HCl	-0.05
Hydrogen Sulfide	H <sub>2</sub> S	-5.0
Nitric Oxide	NO	None
Sulfur Dioxide	SO <sub>2</sub>	-0.8
Phosgene	COCl <sub>2</sub>	None
Chlorine	Cl <sub>2</sub>	1.0

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for a nitrogen dioxide sensor. Note that this response time can become significantly slower at temperatures below -20°C.



## HYDROGEN SENSOR

Hydrogen gas sensor is an electrochemical device used for the measurement of H<sub>2</sub> in ambient air. It is designed to be used in conjunction with Zysense Multi-Gas Analyzer MGA -800. .

H<sub>2</sub> sensors operate by generating a small electrical current proportional to the volumetric concentration of hydrogen gas in the surrounding air. The current is the result of the oxidation of hydrogen on the surface of a catalytic electrode, with a resulting signal that is linear with respect to hydrogen concentration. Hydrogen sensors are 3-electrode sensors and require the presence of at least 5% oxygen to operate properly. They are not suitable for use in oxygen free environments.




Hydrogen sensors respond to very few other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any potential interferences are present and might present interference issues. Note that cross-sensitivity data is approximate and based on exposures under 100 PPM. In some cases, response to other gases may not be stable or may be transient.


### Hydrogen Sensor Cross Sensitivity Data

Gas	Symbol	Response to 1 PPM
Carbon Monoxide	CO	0.1
Acetylene	C <sub>2</sub> H <sub>2</sub>	0.1
Ethylene	C <sub>2</sub> H <sub>4</sub>	0.1
Ethanol (alcohol)	C <sub>2</sub> H <sub>6</sub> O	None
Nitrogen Dioxide	NO <sub>2</sub>	None
Sulfur Dioxide	SO <sub>2</sub>	None
Hydrogen Sulfide	H <sub>2</sub> S	None
Methyl Mercaptan	CH <sub>3</sub> SH	None
Chlorine	Cl <sub>2</sub>	None
Hydrogen Cyanide	HCN	None
Hydrogen Peroxide	H <sub>2</sub> O <sub>2</sub>	None
Ammonia	NH <sub>3</sub>	None
Hydrogen Chloride	HCl	None
Carbon Dioxide	CO <sub>2</sub>	None
Nitric Oxide	NO	None
Hydrogen Fluoride	HF	None

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span of changing temperature



Shown below is a typical response time graph for an H<sub>2</sub> sensor. Note that this response time can become significantly slower at temperatures below -25°C.



## PHOSPHINE SENSOR - PH<sub>3</sub>

Zysense Phosphine sensor is an electrochemical device used for the measurement of very low levels of phosphine gas in ambient air. It is primarily used for parts-per-billion level leak detection. The sensor is actually a “hydride” sensor and is used for other hydride gases such as arsine, diborane, and hydrogen selenide. The cross-sensitivity table on page 2 lists the relative response of this sensor to other hydride gases. Zysense Phosphine sensor is designed to be used in conjunction with Zysense Multi-Gas Analyzer MGA-800.

The sensor is a 3-electrode electrochemical cell and operates by generating a small electrical current proportional to the partial pressure of hydride gas in the surrounding air. The current results from the direct oxidation of phosphine gas on the measuring electrode as shown below.




H10-32 PPB Phosphine sensors exhibit response to few other gases and vapors. When applying this sensor to specific applications, the user should verify whether vapors other than the target vapor are present. Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.

### Phosphine Sensor Cross Sensitivity Data


Gas	Symbol	Response to 1 PPM
Phosphine	PH <sub>3</sub>	1.0
Arsine	AsH <sub>3</sub>	1.0
Diborane	B <sub>2</sub> H <sub>6</sub>	1.5
Hydrogen Selenide	H <sub>2</sub> Se	0.7
Germane	GeH <sub>4</sub>	0.5
Silane	SiH <sub>4</sub>	0.005
Hydrogen Cyanide	HCN	0.05
Acetylene	C <sub>2</sub> H <sub>2</sub>	0.0001
Chlorine (Halogens)	Cl <sub>2</sub>	-0.01
Nitrogen Dioxide	NO <sub>2</sub>	-0.001
Sulfur Dioxide	SO <sub>2</sub>	None
Ethylene	C <sub>2</sub> H <sub>4</sub>	None
Carbon Monoxide	CO	None
Nitric Oxide	NO	None
Alcohol	C <sub>2</sub> H <sub>6</sub> O	None
Hydrogen	H <sub>2</sub>	None



Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for a Phosphine. Response to other hydride gases may be somewhat different. Note that response times can become significantly slower at temperatures below -20°C.



## HYDROGEN SULFIDE SENSOR (H2S)

Model H10-24 Hydrogen Sulfide sensor is an electrochemical device used for the detection of H<sub>2</sub>S gas leaks in ambient air. It is designed to be used in conjunction with Zysense MGA 800 Multi-Gas Analyzer. . Hydrogen Sulfide sensors contain internal electronics and memory that control sensor bias and store calibration data, calibration history, and limited data log. H<sub>2</sub>S sensors operate by generating a small electrical current proportional to the partial pressure of H<sub>2</sub>S gas in the surrounding air. The current is the result of the direct oxidation of sulfide to sulfate as shown in the equation below. H<sub>2</sub>S sensors are 3-electrode sensors and require oxygen to function. They may not be used in applications where O<sub>2</sub> levels are below 5%.




Hydrogen Sulfide sensors exhibit response to certain other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any of these potential interferences are present and might present interference issues.

Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.


### Hydrogen Sulfide Sensor Cross Sensitivity Data

Gas	Symbol	Response to 1 PPM
Ammonia	NH <sub>3</sub>	-0.005
Hydrogen Chloride	HCl	-0.005
Carbon Monoxide	CO	0.01
Carbon Dioxide	CO <sub>2</sub>	None
Nitric Oxide	NO	0.40
Ozone	O <sub>3</sub>	-0.05
Chlorine	Cl <sub>2</sub>	-0.05
Hydrogen	H <sub>2</sub>	0.0003
Nitrogen Dioxide	NO <sub>2</sub>	0.01
Sulfur Dioxide	SO <sub>2</sub>	0.08
Methyl Mercaptan	CH <sub>3</sub> SH	0.3
Hydrides	PH <sub>3</sub>	0.1
Ethanol (alcohol)	C <sub>2</sub> H <sub>6</sub> O	0.015

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span ( $\mu\text{A}/\text{PPM}$ ) of changing temperature.



Shown below is a typical response time graph for an H<sub>2</sub>S sensor. Note that this response time can become significantly slower at temperatures below -20°C.



## Ammonia Sensor ( NH3)

Ammonia sensor is an electrochemical device used for the detection of NH<sub>3</sub> gas leaks in ambient air. It is designed to be used in conjunction with Zysense multi-gas analyzer MGA -800.

Ammonia sensors operate by generating a small electrical current proportional to the partial pressure of NH<sub>3</sub> gas in the surrounding gas. The current is the result of the electrochemical oxidation of NH<sub>3</sub> on a catalytic electrode as shown below. Ammonia sensors are 3-electrode sensors and require oxygen levels above 5% to operate properly. They may not be used in oxygen free applications.




Ammonia sensors exhibit response to certain other gases. When applying this sensor to specific applications, it is good practice to verify whether or not any of these potential interferences are present and might present interference issues. Note that cross-sensitivity data is approximate. In some cases, response to other gases may not be stable or may be transient.

### Ammonia Sensor Cross Sensitivity Data:

Gas	Symbol	Response to 1 PPM
Dimethylamine	(CH <sub>3</sub> ) <sub>2</sub> NH	0.71
Hydrogen Chloride	HCl	-0.5
Carbon Monoxide	CO	0.01
Carbon Dioxide	CO <sub>2</sub>	None
Nitric Oxide	NO	0.08
Hydrogen	H <sub>2</sub>	0.002
Hydrogen Sulfide	H <sub>2</sub> S	1.0
Nitrogen Dioxide	NO <sub>2</sub>	-0.08
Sulfur Dioxide	SO <sub>2</sub>	-0.1
Methyl Mercaptan	CH <sub>3</sub> SH	0.3
Chlorine	Cl <sub>2</sub>	-0.5
Hydrogen Cyanide	HCN	0.1
Ethanol (alcohol)	C <sub>2</sub> H <sub>6</sub> O	0.015

Electrochemical sensors exhibit a response that is temperature dependent to a limited extent. Although the effect of temperature is not large, it is useful to be aware of the effect. Shown below is a graph showing the effect on span (uA/PPM) of changing temperature.



Shown below is a typical response time graph for an ammonia sensor. Note that this response time can become significantly slower at temperatures below -20°C.

